

Micro-Introducer Kits

IT ALL STARTS HERE.

Successful procedures begin with quality insertion components.

Micro-introducer kits from AngioDynamics include the finest tools available for .018" guidewire placement during vascular procedures.



> THE NEEDLE

Choose from A or B bevel, or super sharp needle options. An echogenic tip is available for higher visibility under ultrasound guidance.



> THE MICRO-INTRODUCER

AngioDynamics offers micro-introducers in 4 and 5 French, with both sizes available in standard or stiffened options. AngioDynamics micro-introducers feature a quarter-twist cam connection that saves time by verifying a secure lock, while seamless dilator-to-sheath transitions promote optimum performance. Choose the stiffened dilator to assist smooth entry into difficult vessels.




> THE GUIDEWIRE

The most important component of the micro-introducer kit is the guidewire. AngioDynamics has a range of atraumatic configurations to meet your specific needs, including the following options:

- Stainless steel wire and tip
- Stainless steel wire with tungsten tip for increased radiopacity
- Nitinol wire with tungsten tip for increased radiopacity plus high kink resistance

The .018" wire tips are atraumatic, designed to deflect off the vessel wall.



ANGIODYNAMICS®

MICRO-INTRODUCER KITS

Description	Catalog #			
	4F	5F	4F stiffened	5F stiffened
45cm stainless steel wire with stainless tip - 7cm/A bevel	□ 06597001	□ 06597002	□ 06597036	□ 06597037
45cm stainless steel wire with stainless tip - 7cm/B bevel	□ 06597003	□ 06597004	□ 06597038	□ 06597039
45cm stainless steel wire with stainless tip - 7cm/super sharp needle	□ 06597005	□ 06597006	□ 06597040	□ 06597041
60cm stainless steel wire with stainless tip - 7cm/B bevel	□ 06597007	□ 06597008	□ 06597042	□ 06597043
60cm stainless steel wire with stainless tip - 7cm/super sharp needle	□ 06597009	□ 06597010	□ 06597044	□ 06597045
45cm stainless steel wire with tungsten tip - 7cm/B bevel	□ 06597011	□ 06597012	□ 06597046	□ 06597047
45cm stainless steel wire with tungsten tip - 7cm/super sharp needle	□ 06597013	□ 06597014	□ 06597048	□ 06597049
45cm nitinol wire with tungsten tip - 7cm/B bevel	□ 06597015	□ 06597016	□ 06597050	□ 06597051
45cm nitinol wire with tungsten tip - 7cm/super sharp needle	□ 06597017	□ 06597018	□ 06597052	□ 06597053
45cm nitinol wire with tungsten tip - 7cm/A bevel	□ 06597019	□ 06597020	□ 06597054	□ 06597055
45cm nitinol wire with tungsten tip - 4cm/B bevel	□ 06597021	–	□ 06597056	–
60cm stainless steel wire with stainless tip - 7cm/B bevel*	–	□ 06597022	–	□ 06597057
60cm stainless steel wire with stainless tip - 7cm/super sharp needle*	–	□ 06597023	–	□ 06597058
45cm stainless steel wire with stainless tip - 7cm/B bevel*	□ 06597024	□ 06597025	□ 06597059	□ 06597060
45cm stainless steel wire with stainless tip - 7cm/super sharp needle*	□ 06597026	□ 06597027	□ 06597061	□ 06597062
45cm stainless steel wire with tungsten tip - 7cm/B bevel*	□ 06597028	□ 06597029	□ 06597063	□ 06597064
45cm stainless steel wire with tungsten tip - 7cm/super sharp needle*	□ 06597030	□ 06597031	□ 06597065	□ 06597066
45cm nitinol wire with tungsten tip - 7cm/B bevel*	□ 06597032	□ 06597033	□ 06597067	□ 06597068
45cm nitinol wire with tungsten tip - 7cm/super sharp needle*	□ 06597034	□ 06597035	□ 06597069	□ 06597070

Each kit includes: one 21 gauge entry needle; one .018" guidewire; one 10cm micro-introducer

* non-echogenic

MICRO-INTRODUCER SETS

Description	Catalog #			
	4F	5F	4F stiffened	5F stiffened
Introducer only	□ 06597101	□ 06597102	□ 06597103	□ 06597104

INDIVIDUAL COMPONENTS

Description	Catalog #
21 gauge x 7cm A bevel introducer needle with echogenic tip	□ 06506801
21 gauge x 7cm B bevel introducer needle with echogenic tip	□ 06506802
21 gauge x 4cm pediatric B bevel introducer needle with echogenic tip	□ 06506803
21 gauge x 7cm B bevel introducer needle with non-echogenic tip	□ 06506804
21 gauge x 4cm pediatric B bevel introducer needle with non-echogenic tip	□ 06506805
21 gauge x 7cm super sharp needle with echogenic tip	□ 06508601
21 gauge x 7cm super sharp needle with non-echogenic tip	□ 06508602
.018" x 45cm stainless steel guidewire with stainless tip	□ 06507201
.018" x 60cm stainless steel guidewire with stainless tip	□ 06507204
.018" x 45cm nitinol guidewire with tungsten tip	□ 06597201
.018" x 45cm stainless steel guidewire with tungsten tip	□ 06597202

>> 10 units per box

ANGIODYNAMICS®

USA > 603 Queensbury Avenue, Queensbury, NY 12804 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360

International > AngioDynamics UK Ltd > Building 2000, Beach Drive, IQ Cambridge, Waterbeach, Cambridge, CB25 9TE > United Kingdom

tel: +44 1223 729 300 > fax: +44 1223 729 349

www.angiodynamics.com

AngioDynamics is a registered trademark of AngioDynamics, Inc.
CAUTION: Federal Law (U.S.) restricts the sale of this device by or on the order of a physician.



Proudly Distributed by
MEDICAL SPECIALTIES AUSTRALIA PTY LIMITED
54 GIBBES STREET, CHATSWOOD NSW 2067
PHONE 02 9417 7955 FAX 02 9417 5779
E-Mail: custserv@msa.com.au
www.msa.com.au

© 2009 AngioDynamics, Inc.

ADC1007-06/2011