



**Total Abscession<sup>®</sup>**  
DRAINAGE CATHETERS


**ANGIODYNAMICS<sup>®</sup>**

# Total Abscession®

## DRAINAGE CATHETERS

### Security

- The Vault® locking mechanism securely fixes the pigtail and prevents tampering or accidental removal
- “The Click” provides audible confirmation that the lock is engaged




### Physician Performance

- Kink resistant shaft recovers rapidly, even if severely bent, knotted or twisted
- Shaft diameter equals the inner diameter of the catheter hub for freer flow and clog-free drainage
- Thermal molded tip for less buckling and kinking upon insertion
- Lubricious Blue Silk™ finish provides easier insertion and pushability
- Aspiration can be accomplished in both the locked and unlocked positions

### Patient Comfort

- Low-profile hub and soft, pliable strain relief provide optimal patient comfort regardless of anatomical access point
- Soft shaft design for a more comfortable patient fit



## Three Different Styles of Drain to Meet Specific Clinical Drainage Needs:



### GENERAL DRAIN

>> Recommended for percutaneous drainage of fluids.



AngioDynamics' Total Abscession General drain features six side holes and is available with or without a radiopaque marker.

### BILIARY DRAIN

>> Recommended for percutaneous transhepatic biliary drainage.



AngioDynamics' Total Abscession Biliary drain features additional side holes on the shaft arranged in a spiral fashion for optimal bile duct drainage. It also has a radiopaque marker to better identify the placement of the catheter on imaging equipment.

### NEPHROSTOMY DRAIN

>> Recommended for percutaneous drainage of the kidneys.



AngioDynamics' Total Abscession Nephrostomy drain features six side holes, and has a tighter curl at the distal end of the catheter to ensure the drain tip lies completely inside of the pigtail avoiding irritation of the renal pelvis.

## GENERAL DRAIN

Description	Catalog #
6.5F x 15cm Total Abscession® Drainage Catheter	□ 14000801
6.5F x 30cm Total Abscession Drainage Catheter	□ 14000802
8F x 30cm Total Abscession Drainage Catheter	□ 14000803
10F x 30cm Total Abscession Drainage Catheter	□ 14000804
2F x 30cm Total Abscession Drainage Catheter	□ 14000805
14F x 30cm Total Abscession Drainage Catheter	□ 14000806
8F x 45cm Total Abscession Drainage Catheter	□ 14000807
10F x 45cm Total Abscession Drainage Catheter	□ 14000808
12F x 45cm Total Abscession Drainage Catheter	□ 14000809
14F x 45cm Total Abscession Drainage Catheter	□ 14000810
8F x 30cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000823
10F x 30cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000824
12F x 30cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000825
14F x 30cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000826
8F x 45cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000827
10F x 45cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000828
12F x 45cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000829
14F x 45cm Total Abscession Drainage Catheter with Radiopaque Marker	□ 14000830

## BILIARY DRAIN

Description	Catalog #
8F x 40cm Total Abscession Biliary Drainage Catheter	□ 14000901
10F x 40cm Total Abscession Biliary Drainage Catheter	□ 14000902
12F x 40cm Total Abscession Biliary Drainage Catheter	□ 14000903
14F x 40cm Total Abscession Biliary Drainage Catheter	□ 14000904

## NEPHROSTOMY DRAIN

Description	Catalog #
6.5F x 30cm Total Abscession Nephrostomy Drainage Catheter	□ 14001001
8F x 30cm Total Abscession Nephrostomy Drainage Catheter	□ 14001002
10F x 30cm Total Abscession Nephrostomy Drainage Catheter	□ 14001003
12F x 30cm Total Abscession Nephrostomy Drainage Catheter	□ 14001004
14F x 30cm Total Abscession Nephrostomy Drainage Catheter	□ 14001005
8F x 45cm Total Abscession Nephrostomy Drainage Catheter	□ 14001006
10F x 45cm Total Abscession Nephrostomy Drainage Catheter	□ 14001007
12F x 45cm Total Abscession Nephrostomy Drainage Catheter	□ 14001008
14F x 45cm Total Abscession Nephrostomy Drainage Catheter	□ 14001009

>> Quantity of 5 per box.

**ANGIODYNAMICS®**

USA > 603 Queensbury Avenue, Queensbury, NY 12804 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360  
 International > AngioDynamics UK Ltd > Building 2000, Beach Drive, IQ Cambridge, Waterbeach, Cambridge, CB25 9TE > United Kingdom  
 tel: +44 1223 729 300 > fax: +44 1223 729 349

[www.angiodynamics.com](http://www.angiodynamics.com)

Abscession, Vault and AngioDynamics are registered trademarks of AngioDynamics, Inc.  
 Blue Silk is a trademark of AngioDynamics.  
 Covered under U.S. Patent 7,217,256  
 CAUTION: Federal Law (U.S.) restricts the sale of this device by or on the order of a physician.

PHYSICIAN'S SIGNATURE



Proudly distributed by:  
**MEDICAL SPECIALTIES AUSTRALIA PTY LIMITED**  
 54 GIBBES STREET, CHATSWOOD NSW 2067  
 PHONE 02 9417 7955 FAX 02 9417 5779  
 E-Mail: [custserv@msa.com.au](mailto:custserv@msa.com.au)  
[www.msa.com.au](http://www.msa.com.au)

© 2009 AngioDynamics, Inc. MLC 246 Rev B 07/09

ADC1014-06/2011