


Plastic Surgery
TiLOOP® Bra

The worldwide unique, extra-light
titanized mesh for reconstructive
breast surgery.




www.pfmmedical.com

titanized polypropylene



TiLOOP® Bra

- the worldwide unique, extra light titanized mesh for reconstructive breast surgery.



Benefits

Titanization

- ▶ Highly biocompatible material
- ▶ Excellent ingrowth
- ▶ Breast tissue stays soft
- ▶ Women keep their natural breast feeling

Details

Characteristics

- ▶ Made of titanized polypropylene
- ▶ Extra - light with just 16 g/m² of weight
- ▶ 3 different sizes
- ▶ No influence on MRI. X-ray and ultrasound

Indications

- ▶ Primary breast reconstruction
- ▶ Secondary breast reconstruction
- ▶ Implant exchange
- ▶ Mastopexy

Ordering information

Name	Size	REF	Package Unit
TiLOOP® Bra extralight	small	6000636	1
TiLOOP® Bra extralight	medium	6000637	1
TiLOOP® Bra extralight	large	6000638	1



MEDICAL
SPECIALTIES
AUSTRALIA

54 Gibbes St, Chatswood 2067 NSW, Australia
PO Box 764, Willoughby 2068 NSW, Australia

✉ custserv@msa.com.au ☎ 02 9417 7955
🌐 www.msa.com.au ☎ 02 9417 5779



pfm medical ag

Wankelstraße 60 50996 Köln, Germany

Phone +49 (0) 2236 96 41-0

Fax +49 (0) 2236 96 41- 744

PUE1012-06/2014