


Soft-Vu[®]





ANGIOGRAPHIC CATHETERS











ANGIODYNAMICS[®]

**Cerebral
SELECTIVE
CATHETERS**

Soft Radiopaque Tip





	 Sidewinder ST 1	 Sidewinder ST 2	 Straight Art.	 Straight Art.				
	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.
5 per Box	4F / 65 cm / .038"	□ 10714030	4F / 65 cm / .038"	□ 10714031	*4F / 65 cm / .035"	□ 10731402	4F / 65 cm / .035"	□ 10731401
High-Torque Shaft	4F / 100 cm / .038"	□ 10714041	4F / 100 cm / .038"	□ 10714042	*4F / 90 cm / .035"	□ 10714015	4F / 90 cm / .035"	□ 10731403
*1.5 cm tip					*5F / 65 cm / .035"	□ 10714022	5F / 65 cm / .035"	□ 10731502
**2 cm tip					**5F / 25 cm / .035"	□ 10731501	5F / 25 cm / .035"	□ 10731509
							5F / 40 cm / .035"	□ 10731510
Braided Ultra-Torque Shaft	4F / 100 cm / .035"	□ 10707605	4F / 100 cm / .035"	□ 10707606				
	5F / 100 cm / .038"	□ 10709009	5F / 100 cm / .035"	□ 10709011				
			5F / 100 cm / .038"	□ 10709010				
					High-Torque Shaft		6F / 100 cm / .054"	□ 10731602
							6F / 110 cm / .054"	□ 10731601
					Braided Ultra-Torque Shaft		4F / 130 cm / .035"	□ 10731503





	 Bentson 1	 Bentson 2	 Bentson 3	 Cherne	
High-Torque Shaft	4F / 100 cm / .035"	□ 10719001	4F / 100 cm / .035"	□ 10719002	
	5F / 100 cm / .035"	□ 10720201	4F / 65 cm / .038"	□ 10714028	
	5F / 100 cm / .038"	□ 10720204	4F / 100 cm / .038"	□ 10714043	
			5F / 100 cm / .035"	□ 10720202	
			5F / 100 cm / .038"	□ 10720205	
Braided Ultra-Torque Shaft	4F / 100 cm / .035"	□ 10708301	4F / 100 cm / .035"	□ 10708302	
	5F / 100 cm / .035"	□ 10709701	5F / 100 cm / .035"	□ 10709702	
	5F / 100 cm / .038"	□ 10709704	5F / 100 cm / .038"	□ 10709705	
	6F / 100 cm / .038"	□ 10711104	6F / 100 cm / .035"	□ 10711102	
			6F / 100 cm / .038"	□ 10711105	
				5F / 69 cm / .035"	□ 10735201





	 Head Hunter 1	 Head Hunter 3	 JB 1	 JB 2
High-Torque Shaft	4F / 100 cm / .035"	□ 10718301	4F / 100 cm / .035"	□ 10718302
*1 cm tip	4F / 65 cm / .038"	□ 10714027	5F / 100 cm / .035"	□ 10719502
**3 cm tip	4F / 100 cm / .038"	□ 10714040	5F / 100 cm / .038"	□ 10719504
***8.5 cm tip	5F / 100 cm / .035"	□ 10719501	6F / 100 cm / .035"	□ 10720702
	5F / 100 cm / .038"	□ 10719503	5F / 130 cm / .038"	□ 10734208
	5F / 140 cm / .038"	□ 10719505		
	6F / 100 cm / .035"	□ 10720701		
	6F / 100 cm / .038"	□ 10720703		
Braided Ultra-Torque Shaft	4F / 100 cm / .035"	□ 10707501	4F / 100 cm / .035"	□ 10707502
**3 cm tip	5F / 100 cm / .035"	□ 10708901	5F / 100 cm / .035"	□ 10708902
****5 cm tip	5F / 65 cm / .038"	□ 10708905	5F / 100 cm / .038"	□ 10708904
	5F / 100 cm / .038"	□ 10708903		
	4F / 130 cm / .035"	□ 10708906	**4F / 65 cm / .035"	□ 10734001
			**4F / 100 cm / .038"	□ 10714049
			**5F / 100 cm / .035"	□ 10734201
			*5F / 100 cm / .035"	□ 10734204
			**5F / 100 cm / .035"	□ 10734205
			****5F / 110cm / .035"	□ 10734203
			*5F / 130 cm / .038"	□ 10734207
			**5F / 100 cm / .035"	□ 10734103
			**5F / 65 cm / .038"	□ 10734101
			**5F / 100 cm / .038"	□ 10734102

Cerebral
**SELECTIVE
CATHETERS**

Soft Radiopaque Tip

	 JB 3		 Kerber		 Newton 1		 Newton 2	
	Description	Catalog	Description	Catalog	Description	Catalog No.	Description	Catalog No.
5 per Box								
High-Torque Shaft	*4F / 65 cm / .038"	□ 10714037	5F / 100 cm / .035"	□ 10734801	4F / 100 cm / .035"	□ 10719101	4F / 100 cm / .035"	□ 10719102
*3 cm tip	*4F / 100 cm / .038"	□ 10714051			5F / 100 cm / .035"	□ 10720301	5F / 100 cm / .035"	□ 10720302
**4 cm tip	**5F / 100 cm / .035"	□ 10734206			5F / 100 cm / .038"	□ 10720304	5F / 100 cm / .038"	□ 10720305
Braided Ultra-Torque Shaft			4F / 100 cm / .035"	□ 10734601	4F / 100 cm / .035"	□ 10708401	4F / 100 cm / .035"	□ 10708402
			5F / 100 cm / .038"	□ 10734701	5F / 100 cm / .035"	□ 10709801	5F / 100 cm / .035"	□ 10709802
					5F / 100 cm / .038"	□ 10709804	5F / 100 cm / .038"	□ 10709805



	 Newton 3		 Newton 5		 Mani		 Osborne 1	
	Description	Catalog	Description	Catalog	Description	Catalog No.	Description	Catalog No.
High-Torque Shaft	4F / 100 cm / .035"	□ 10719103	4F / 100 cm / .035"	□ 10719104	4F / 100 cm / .035"	□ 10718201	4F / 100 cm / .035"	□ 10732902
	5F / 100 cm / .035"	□ 10720303	4F / 65 cm / .038"	□ 10714032	4F / 65 cm / .038"	□ 10714029		
	5F / 100 cm / .038"	□ 10720306	4F / 100 cm / .038"	□ 10714046	4F / 100 cm / .038"	□ 10714044		
			5F / 100 cm / .035"	□ 10720307	5F / 100 cm / .035"	□ 10719401		
			5F / 140 cm / .038"	□ 10720308	5F / 100 cm / .038"	□ 10719402		
Braided Ultra-Torque Shaft	4F / 100 cm / .035"	□ 10708403			4F / 100 cm / .035"	□ 10707401	4F / 100 cm / .035"	□ 10734501
	5F / 100 cm / .035"	□ 10709803			5F / 100 cm / .035"	□ 10708801		
	5F / 100 cm / .038"	□ 10709806			5F / 100 cm / .038"	□ 10708802		

	 Osborne 2		 Sidewinder 1		 Sidewinder 2		 Sidewinder 3	
	Description	Catalog	Description	Catalog	Description	Catalog No.	Description	Catalog No.
High-Torque Shaft	5F / 100 cm / .035"	□ 10734401	4F / 65 cm / .035"	□ 10718401	4F / 100 cm / .035"	□ 10718403	4F / 100 cm / .035"	□ 10718404
	4F / 100 cm / .035"	□ 10732901	4F / 100 cm / .035"	□ 10718402	5F / 100 cm / .035"	□ 10719603	5F / 100 cm / .035"	□ 10719604
			5F / 65 cm / .035"	□ 10719601	5F / 100 cm / .038"	□ 10719607	5F / 100 cm / .038"	□ 10719608
			5F / 100 cm / .035"	□ 10719602	6F / 100 cm / .035"	□ 10710403	6F / 100 cm / .035"	□ 10710404
			5F / 65 cm / .038"	□ 10719605	6F / 100 cm / .038"	□ 10710407		
			5F / 100 cm / .038"	□ 10719606				
			4F / 65 cm / .038"	□ 10719609				
			6F / 100 cm / .035"	□ 10720802				
			6F / 65 cm / .035"	□ 10720801				
Braided Ultra-Torque Shaft			4F / 65 cm / .035"	□ 10707601	4F / 100 cm / .035"	□ 10707603	4F / 100 cm / .035"	□ 10707604
			4F / 100 cm / .035"	□ 10707602	5F / 100 cm / .035"	□ 10709003	5F / 100 cm / .035"	□ 10709004
			5F / 65 cm / .035"	□ 10709001	5F / 100 cm / .038"	□ 10709007	5F / 100 cm / .038"	□ 10709008
			5F / 100 cm / .035"	□ 10709002	4F / 100 cm / .035"	□ 10707606		
			5F / 65 cm / .038"	□ 10709005	5F / 100 cm / .038"	□ 10709010		
			5F / 100 cm / .038"	□ 10709006	5F / 100 cm / .035"	□ 10709013		
			6F / 100 cm / .038"	□ 10710406				
			6F / 65 cm / .035"	□ 10710401				
			6F / 65 cm / .038"	□ 10710405				

**Angled
SELECTIVE
CATHETERS**

Soft Radiopaque Tip

5 per Box




High-Torque Shaft

Description	Catalog No.	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.
4F / 65 cm / .038"	□ 10714033	4F / 65 cm / .038"	□ 10714034	4F / 65 cm / .035"	□ 10733801	4F / 40 cm / .035"	□ 10734301
4F / 100 cm / .038"	□ 10714047	4F / 100 cm / .038"	□ 10714048			4F / 100 cm / .035"	□ 10734302
5F / 65 cm / .035"	□ 10721805	5F / 65 cm / .035"	□ 10719902			4F / 65 cm / .038"	□ 10714014
5F / 100 cm / .035"	□ 10721801	5F / 100 cm / .035"	□ 10721802			4F / 100 cm / .038"	□ 10714021
5F / 100 cm / .038"	□ 10721803	5F / 100 cm / .038"	□ 10721804			5F / 40 cm / .038"	□ 10732702
5F / 65 cm / .035"	□ 10721805					5F / 65 cm / .038"	□ 10732701
						5F / 100 cm / .038"	□ 10732703
						5F / 25 cm / .035"	□ 10732707

Braided Ultra-Torque Shaft

5F / 100 cm / .035"	□ 10722406	4F / 100 cm / .035"	□ 10723302	4F / 65 cm / .035"	□ 10733901
5F / 100 cm / .038"	□ 10722403	5F / 100 cm / .035"	□ 10722402		
4F / 100 cm / .035"	□ 10732501	5F / 65 cm / .038"	□ 10722405		
		5F / 100 cm / .038"	□ 10722404		



High-Torque Shaft



*10 cm tip

**w/ 6 side holes

4F / 65 cm / .038"	□ 10714038	4F / 65 cm / .035"	□ 10735101	5F / 65 cm / .035"	□ 10735601	4F / 65 cm / .038"	□ 10714009
4F / 100 cm / .038"	□ 10714052	4F / 65 cm / .038"	□ 10714039	**5F / 100 cm / .035"	□ 10735602	4F / 100 cm / .038"	□ 10714016
6F / 90 cm / .035"	□ 10735701	4F / 100 cm / .038"	□ 10714053			4F / 110 cm / .035"	□ 10723601
		5F / 65 cm / .035"	□ 10735001			5F / 100 cm / .035"	□ 10722101
						5F / 40 cm / .038"	□ 10722105
						*5F / 65 cm / .038"	□ 10722108
						5F / 65 cm / .038"	□ 10722104
						5F / 100 cm / .038"	□ 10722102
						5F / 110 cm / .038"	□ 10722103
						5F / 65 cm / .046"	□ 10722106
						5F / 100 cm / .046"	□ 10722107
						6F / 65 cm / .054"	□ 10731701
						5F / 140 cm / .038"	□ 10722109
						6F / 65 cm / .038"	□ 10731702
						4F / 40 cm / .035"	□ 10723503
						4F / 65 cm / .035"	□ 10723502
						4F / 100 cm / .035"	□ 10723501
						4F / 110 cm / .035"	□ 10723504
						5F / 100 cm / .035"	□ 10722701
						5F / 65 cm / .035"	□ 10722704
						5F / 65 cm / .038"	□ 10722703
						5F / 100 cm / .038"	□ 10722702
						5F / 110 cm / .038"	□ 10722705
						6F / 100 cm / .038"	□ 10711301
						6F / 110 cm / .038"	□ 10711302
						5F / 140 cm / .035"	□ 10722706
						5F / 130 cm / .038"	□ 10722707
						4F / 140 cm / .035"	□ 10722708

Braided Ultra-Torque Shaft

4F / 65 cm / .035"	□ 10735401	5F / 65 cm / .038"	□ 10734901
4F / 100 cm / .035"	□ 10735402		
5F / 50 cm / .035"	□ 10735501		
5F / 80 cm / .035"	□ 10735503		
5F / 100cm / .038"	□ 10735504		



High-Torque Shaft

4F / 100 cm / .038"	□ 10714045	5F / 100 cm / .035"	□ 10733601
4F / 100 cm / .035"	□ 10718101		
5F / 100 cm / .035"	□ 10719301		
6F / 100 cm / .038"	□ 10720502		

Braided Ultra-Torque Shaft

4F / 100 cm / .035"	□ 10707301	5F / 100 cm / .035"	□ 10733502
5F / 100 cm / .035"	□ 10708701	5F / 100 cm / .038"	□ 10733501
5F / 100 cm / .038"	□ 10708702		

Visceral
**SELECTIVE
CATHETERS**

Soft Radiopaque Tip

5 per Box

High-Torque Shaft





*Sos Omni 0 Mini

**Sos Omni 2 Mini

***w/ 2 side holes

Braided Ultra-Torque Shaft

****w/ 4 side holes


	 Sos Omni 0		 Sos Omni 1		 Sos Omni 2		 Sos Omni 3	
	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.
	4F / 80 cm / .038"	□ 10714023	4F / 80 cm / .038"	□ 10714024	4F / 80 cm / .038"	□ 10714025	4F / 80 cm / .038"	□ 10714026
	5F / 80 cm / .035"	□ 10720407	5F / 80 cm / .035"	□ 10720402	5F / 80 cm / .035"	□ 10720403	5F / 80 cm / .035"	□ 10720404
	5F / 80 cm / .038"	□ 10720408	***5F / 80 cm / .035"	□ 10732002	**5F / 80 cm / .035"	□ 10732001	5F / 80 cm / .038"	□ 10720406
	*4F / 65 cm / .035"	□ 10711502	5F / 80 cm / .038"	□ 10720401	5F / 80 cm / .038"	□ 10720405		
	4F / 80 cm / .035"	□ 10711404	4F / 80 cm / .035"	□ 10711401	4F / 80 cm / .038"	□ 10711402	4F / 80 cm / .035"	□ 10711403
	5F / 80 cm / .035"	□ 10712004	5F / 80 cm / .035"	□ 10712001	****4F / 80 cm / .035"	□ 10711405	5F / 80 cm / .035"	□ 10712003
	5F / 80 cm / .038"	□ 10712008	5F / 80 cm / .038"	□ 10712005	5F / 80 cm / .035"	□ 10712002	5F / 80 cm / .038"	□ 10712007
					5F / 80 cm / .038"	□ 10712006		

High-Torque Shaft

*w/ 2 side holes





Braided Ultra-Torque Shaft

*w/ 2 side holes

	 Hook 1		 Cobra 1		 Cobra 2		 Cobra 3	
	4F / 65 cm / .035"	□ 10718701	4F / 65 cm / .035"	□ 10718501	4F / 65 cm / .035"	□ 10718502	4F / 65 cm / .035"	□ 10718503
	5F / 65 cm / .035"	□ 10719901	4F / 90 cm / .035"	□ 10718504	4F / 90 cm / .035"	□ 10718505	4F / 90 cm / .035"	□ 10718506
	5F / 65 cm / .038"	□ 10719903	4F / 65 cm / .038"	□ 10714010	4F / 65 cm / .038"	□ 10714011	4F / 65 cm / .038"	□ 10714012
			4F / 100 cm / .038"	□ 10714017	4F / 100 cm / .038"	□ 10714018	4F / 100 cm / .038"	□ 10714019
			5F / 65 cm / .035"	□ 10719701	5F / 65 cm / .035"	□ 10719702	5F / 90 cm / .035"	□ 10719706
			5F / 90 cm / .035"	□ 10719704	5F / 90 cm / .035"	□ 10719705	5F / 65 cm / .038"	□ 10719709
			5F / 65 cm / .038"	□ 10719707	5F / 65 cm / .038"	□ 10719708	5F / 90 cm / .038"	□ 10719712
			5F / 90 cm / .038"	□ 10719710	5F / 90 cm / .038"	□ 10719711	5F / 65 cm / .038"	□ 10731803
			*5F / 65 cm / .038"	□ 10731801	5F / 65 cm / .038"	□ 10731802	5F / 65 cm / .035"	□ 10731804
					5F / 65 cm / .035"	□ 10731805	5F / 65 cm / .035"	□ 10719703
					6F / 90 cm / .035"	□ 10720905	6F / 65 cm / .038"	□ 10720909
					6F / 65 cm / .038"	□ 10720908		
	4F / 65 cm / .035"	□ 10707901	4F / 65 cm / .035"	□ 10707701	4F / 65 cm / .035"	□ 10707702	4F / 65 cm / .035"	□ 10707703
	5F / 65 cm / .035"	□ 10709301	4F / 90 cm / .035"	□ 10707704	4F / 65 cm / .035"	□ 10707708	4F / 90 cm / .035"	□ 10707706
	5F / 65 cm / .038"	□ 10709303	*4F / 65 cm / .035"	□ 10707707	4F / 90 cm / .035"	□ 10707705	5F / 65 cm / .035"	□ 10709103
			5F / 65 cm / .035"	□ 10709101	5F / 65 cm / .035"	□ 10709102	5F / 90 cm / .035"	□ 10709106
			5F / 90 cm / .038"	□ 10709104	5F / 90 cm / .035"	□ 10709105	5F / 65 cm / .038"	□ 10709109
			5F / 65 cm / .038"	□ 10709107	5F / 65 cm / .038"	□ 10709108	5F / 90 cm / .038"	□ 10709112
			5F / 90 cm / .038"	□ 10709110	5F / 90 cm / .038"	□ 10709111		
			6F / 90 cm / .038"	□ 10710510	5F / 80 cm / .038"	□ 10709113		
					6F / 90 cm / .038"	□ 10710511		

High-Torque Shaft

Braided Ultra-Torque Shaft

	 Hook 2		 Shepherd Hook 1		 Shepherd Hook 2		 Renal	
	4F / 65 cm / .035"	□ 10718702	4F / 65 cm / .035"	□ 10718601	4F / 65 cm / .035"	□ 10718602	4F / 65 cm / .035"	□ 10718801
	5F / 65 cm / .035"	□ 10719902	5F / 65 cm / .035"	□ 10719801	5F / 65 cm / .035"	□ 10719802	5F / 65 cm / .035"	□ 10720001
	5F / 65 cm / .038"	□ 10719904	5F / 65 cm / .038"	□ 10719803	5F / 65 cm / .038"	□ 10719804	5F / 65 cm / .038"	□ 10720002
	6F / 65 cm / .035"	□ 10721102	6F / 65 cm / .038"	□ 10721101	6F / 65 cm / .035"	□ 10721002		
	6F / 65 cm / .038"	□ 10721104	6F / 65 cm / .035"	□ 10721001				
	4F / 65 cm / .035"	□ 10707902	4F / 65 cm / .035"	□ 10707801	4F / 65 cm / .035"	□ 10707802	4F / 65 cm / .035"	□ 10708001
	5F / 65 cm / .035"	□ 10709302	5F / 65 cm / .035"	□ 10709201	5F / 65 cm / .035"	□ 10709202	5F / 65 cm / .035"	□ 10709401
	5F / 65 cm / .038"	□ 10709304	5F / 65 cm / .038"	□ 10709203	5F / 65 cm / .038"	□ 10709204	5F / 65 cm / .038"	□ 10709402
			6F / 65 cm / .038"	□ 10710603				

Visceral SELECTIVE CATHETERS

Soft Radiopaque Tip



5 per Box	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.
Braided Ultra-Torque Shaft	5F / 40 cm / .035"	□ 10723042	4F / 90 cm / .035"	□ 10733701	4F / 65 cm / .035"	□ 10733001	5F / 65 cm / .038"	□ 10733102
	5F / 65 cm / .035"	□ 10723018	5F / 80 cm / .038"	□ 10732101	5F / 65 cm / .035"	□ 10733101		
	5F / 25 cm / .035"	□ 10723074	4F / 100 cm / .038"	□ 10714044				



High-Torque Shaft	Description	Catalog No.	Description	Catalog No.
	5F / 65 cm / .038"	□ 10733401		
Braided Ultra-Torque Shaft	4F / 65 cm / .035"	□ 10732801	4F / 65 cm / .035"	□ 10733201
			5F / 65 cm / .035"	□ 10733301
			5F / 40 cm / .035"	□ 10733402
			5F / 25 cm / .035"	□ 10733403

FLUSH CATHETERS

Soft Radiopaque Tip

10 per Box



High-Torque Shaft	Description	Catalog No.	Description	Catalog No.	Description	Catalog No.
	4F / 65 cm / .035"	□ 10732301	4F / 65 cm / .035"	□ 10707101	4F / 65 cm / .035"	□ 10707201
*Omni Flush Mini	4F / 90 cm / .035"	□ 10732302	4F / 90 cm / .035"	□ 10707102	4F / 90 cm / .035"	□ 10707202
**Pigtail Mini	4F / 110 cm / .035"	□ 10732303	4F / 100 cm / .035"	□ 10707103	4F / 65 cm / .038"	□ 10714003
	4F / 65 cm / .038"	□ 10714001	4F / 110 cm / .035"	□ 10707104	4F / 90 cm / .038"	□ 10714006
	4F / 90 cm / .038"	□ 10714004	4F / 65 cm / .038"	□ 10714002	5F / 65 cm / .035"	□ 10708601
	4F / 100 cm / .038"	□ 10714007	4F / 90 cm / .038"	□ 10714005	5F / 90 cm / .035"	□ 10708602
	5F / 65 cm / .035"	□ 10732201	4F / 110 cm / .038"	□ 10714008	5F / 65 cm / .038"	□ 10708603
	5F / 90 cm / .035"	□ 10732202	5F / 65 cm / .035"	□ 10708501	5F / 90 cm / .038"	□ 10708604
	5F / 65 cm / .038"	□ 10732203	5F / 90 cm / .035"	□ 10708502	6F / 65 cm / .035"	□ 10710001
	5F / 90 cm / .038"	□ 10732204	5F / 110 cm / .035"	□ 10708503	6F / 90 cm / .035"	□ 10710002
	5F / 100 cm / .038"	□ 10732205	5F / 65 cm / .038"	□ 10708504	6F / 65 cm / .038"	□ 10710003
	*4F / 65 cm / .035"	□ 10732304	5F / 90 cm / .038"	□ 10708505	6F / 90 cm / .038"	□ 10710004
			5F / 110 cm / .038"	□ 10708506		
			5F / 130 cm / .038"	□ 10708507		
			6F / 110 cm / .035"	□ 10709903		
			6F / 110 cm / .038"	□ 10709906		
			**4F / 65 cm / .035"	□ 10707105		

PHYSICIAN'S SIGNATURE
Specialized Catheter Options Available

IMPORTANT RISK INFORMATION

INDICATION FOR USE: AngioDynamics Angiographic Catheters are designed for use where angiographic diagnosis is indicated.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

CONTRAINDICATIONS: Angiographic Catheters should not be used for any purpose other than those indicated in the instructions.

WARNINGS AND PRECAUTIONS: This device is sterilized by ethylene oxide and is intended for single patient use only. Do not reuse or

resterilize. Contents sterile in unopened, undamaged package. Do not use if opened or any sign of product damage is visible.

AngioDynamics Angiographic Catheters should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures. Do not insert catheters directly through synthetic vascular grafts. Insert through a sheath introducer. AngioDynamics Angiographic Catheters are designed for use with specific guidewire diameters. The recommended maximum guidewire diameter is specified on the catheter label.

Optimal guidewire size and judicious use are recommended. Please see package insert for complete list of warnings and precautions.

POTENTIAL COMPLICATIONS: The following adverse reactions have been reported and are associated with the use of angiographic catheters: Thrombus formation, emboli, arterial wall damage, plaque dislodgment, hematoma, cardiac arrhythmias, myocardial infarction, stroke, and death.

ANGIODYNAMICS®

USA > 603 Queensbury Avenue, Queensbury, NY 12804 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Building 2000, Beach Drive, IQ Cambridge, Waterbeach, Cambridge, CB25 9TE > United Kingdom
tel: +44 1223 729 300 > fax: +44 1223 729 349

www.angiodynamics.com

AngioDynamics, Soft-Vu and SOS Omni are registered trademarks of AngioDynamics, Inc. Omni is a trademark of AngioDynamics, Inc.
U.S Patents: 5,584,821 and 5,651,776.
AngioDynamics®, Inc. All rights reserved.

**MEDICAL SPECIALTIES
AUSTRALIA PTY LIMITED**
54 GIBBES ST, CHATSWOOD NSW 2067
PO BOX 764, WILLOUGHBY, NSW 2068
PHONE 02 9417 7955 FAX 02 9417 5779
E-Mail: custserv@msa.com.au
www.msa.com.au



ADC1044-12/2011

© 2011 AngioDynamics, Inc. All rights reserved. MLC 065 US Rev L 04/11