

# THE POWERWAND™ EDC

## Extended Dwell Catheter

Longer than a PIV...shorter than a midline



### Ideal for:

- Difficult IV access
- US-guided PIVs
- Renal patients
- Vein preservation

A whole new class of VAD



access scientific  
.....

# The new 3Fr, 6cm **THE POWERWAND™ EDC** first in class, best in class

Made of Chronoflex C® – a unique material which allows for the lowest complications and highest completion of therapy of any catheter on the market

## Unrivaled Performance

- ▼ Extended Dwell (up to 29 days)
- ▼ Blood Drawable
- ▼ High Flow (75mL/min)
- ▼ Kink Resistant
- ▼ Power Injectable (8mL/sec)

POWERWAND EDC gives patients what they deserve - the potential for a **one-stick hospitalization**



### Ordering Information:

CODE	DESCRIPTION	QTY/BOX
94222	POWERWAND XL Single Sterile, 3F (6cm)	10

POWERWAND is trademark of Access Scientific, LLC.