



Forty Years of
Swiss Innovation



The Swiss Instruments

Experience the all new Nouvag Infiltration Pump Dispenser DP30 with its high performance pump system. Preset the pump performance and work hands free by using the foot pedal.

Dispenser DP30

From the specialist for sturdy pumps



MEDICAL
SPECIALTIES
AUSTRALIA

54 Gibbes St, Chatswood 2067 NSW, Australia
PO Box 764, Willoughby 2068 NSW, Australia

✉ custserv@msa.com.au ☎ 02 9417 7955
🌐 www.msa.com.au 📠 02 9417 5779

Nouvag AG · St.Gallerstr. 23-25 · CH-9403 Goldach
Phone +41 (0)71 846 66 00 · Fax +41 (0)71 846 66 70
info@nouvag.com · www.nouvag.com

Nouvag GmbH · Schulthaißstr. 15 · D-78462 Konstanz
Phone +49 (0)7531 1290-0 · Fax +49 (0)7531 1290-12
info-de@nouvag.com · www.nouvag.com

Universally usable Infiltration Pump

The Dispenser DP 30 is a universally usable and controllable infiltration pump which is highly suitable for deployment in tumescent applications and in vein treatments.

The percentage infiltration amount is controlled using the rotary potentiometer and the foot pedal.

- Best price/performance ratio
- Universal and flexible applicable
- Easy handling with low maintenance effort
- Time saving compared to time consuming injection of the infiltration quality
- Continuously variable and precise adjustable infiltration performance
- The adipolysis process is decisively improved; subsequent fat suction is significantly easier
- Reliable, approved and stable technology

Technical Data, Dispenser DP 30

Type:	Dispenser DP 30
Mains voltage:	100/115/230 Volt at 50 – 60 Hz
Max. output power:	40 VA
Coolant flow rate:	1 – 17 l/h
Max. pressure:	1.5 bar
Width/depth/height:	260/250/110 mm
Weight, control unit:	2.5 kg

Art.-No. 4180, Control unit

Included:

- Infiltration Tube Set No. 6022, disposable
- Bottle holder Art.-No. 1770 (without bottle)

Optional:

- Foot switch On/Off Art.-No. 1513, or...
- Foot pedal Vario Art.-No. 1501

Infiltration liquids not included!

Ergonomic housing

Lightest and smallest housing of its class.
Simple to clean surfaces.

On/Off Control lamp

A green LED indicates if the device is switched on.

Integrated liquid container

Standard infiltration solution bottle are conveniently suspended on a plug-in mount.

Easy set up of the tube

Wide opening compartment for easy set up of the tube.

Low noise pump drive

The infiltration tube is easy to integrate and the interior pump works at a low noise level.

Accurate pump performance

Reproducible infiltration performance due to scaled regulator.

Front access to the pump

Access to the pump with just one hand is at the press of a button next to the display, hence no inadvertent movement of the device.

